The Pure Truth

法科学DNA分析用消耗品：Eppendorf Forensic DNA Grade according to ISO 18385

»Eppendorf Forensic DNA Grade, for results you can trust.«

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Eppendorf Certificate 

Certificate of Purity – Eppendorf Forensic DNA Grade according to ISO 18385

This package contains a high-quality consumable manufactured under the "Forensic DNA Grade according to ISO 18385" Eppendorf Purity Standard.

The ISO 18385 Forensic DNA Grade consumables are produced in a clean room environment; class 6 (according to VDI 2083) and class 100,000 (according to U.S. Fed. Stand. 209 D).

For this product Eppendorf certifies the following:

Free of detectable

- > Human DNA
- > DNase
- > RNase
- > PCR inhibitors

Quality control and subsequent certification is performed by an independent laboratory accredited according to ISO 17025. Lot-specific certificates are available on request or on the internet at www.eppendorf.com/certificates. The product manual is available at: www.eppendorf.com/manuals

The certification comprises following tests:

Human DNA Contamination Test
A probe-based real-time PCR master mix is prepared for the detection of human DNA. The primers amplify a 62 bp fragment present in more than 1×10^6 copies per human cell. The detection of this fragment is performed with a fluorescently labeled DNA probe. Additionally, primers and DNA probes for detecting an internal positive control (IPC) are also added to the master mix. This master mix is used for running positive control, negative control, and test samples.

Positive control: 10 μ L human DNA (0.5 pg/ μ L) and IPC DNA are added to 15 μ L master mix.
Negative control: 10 μ L human DNA-free H₂O and IPC DNA are added to 15 μ L master mix.
Test sample: 15 consumable samples are rinsed one after another with DNA-free water. As an extraction control, IPC DNA is added to the rinse water prior to DNA extraction. Subsequently, an extraction procedure using the standard protocol of a DNA extraction kit is applied on the rinse water resulting in an eluate of 100 μ L. 10 μ L of this solution are added to 15 μ L master mix.


The emission of a fluorescence signal is detected in samples and controls. For the samples to pass certification, no fluorescence signal of the human DNA probe must be found corresponding to the negative control.

Page 1 of 2 

www.eppendorf.com

> Lot specific certificates available at:
www.eppendorf.com/certificates

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