




# Data Harmony

BioNsisht® cloud – Your data, always available

# Your Data is Simply There

All your bioprocess data available - anytime, with just one click

BioNsight® cloud is a cloud-based software solution for bioprocess monitoring and analysis. It is fully integrated into the Eppendorf bioprocess control software DASware® control 6 and enables you to consolidate your bioprocess data in one central location. This data storage is accessible from everywhere, by everyone you authorize. You can easily merge data from different runs, generated with different bioprocess controllers, to transfer data to knowledge. You and your colleagues are able to collaborate on the same data sets, without the need of being at the same place.



## Securing your data

BioNsight cloud was developed with the highest standards of data security in mind. Established and trusted technologies give you the peace of mind that your precious data is well protected, yet accessible.

- > Hosted and built up on Microsoft® Azure technology
- > Data is encrypted at rest using a 256-bit AES encryption
- > Microsoft SOC2/ISI 27001 certification

## Data export

Your data is not only accessible within BioNsight cloud but can easily be exported.

- > Data export in machine-readable JSON format


# Remote Monitoring

For the peace of mind that your bioprocess is going well

BioNsight cloud empowers you and your team with access to real-time bioprocess data straight from the cloud, whether you're in the lab, in the office or on a business trip. So you always precisely know how your bioprocess runs are performing and can make the right decisions faster, while your process is still running. Or simply sit back with the peace of mind that all is well.

## View all your running processes at a glance

Maintain an instant overview of all your devices to be always informed about the status of your bioprocesses.



## Select one unit for a closer view

### Look back in time

Access all process data since the start of the bioprocess to make informed decisions.

### Get quick overview


Compare current process values and setpoints at a glance.

### Compare with golden batch

Add reference curve from historic process records to quickly spot deviations.

### Select your parameters of interest

Change your selection anytime to display what is relevant to you.



# Integrate your Data

Merge all your datasets at one place to see the entire picture

BioNsight cloud allows you to instantly view all relevant bioprocess data in a vendor-agnostic manner. You can automatically transfer data from internal sensors during a running process, upload offline data from external analyzers, include data from historic runs, and integrate process data from non-Eppendorf bioprocess controllers. This accessibility maximizes the potential of your bioprocess data. Additionally, it makes it easy to grant data access to internal and external partners for discussion, analysis, and troubleshooting.

## Data transfer to the cloud in real-time

- > Connect your existing DASware control 6 installation to BioNsight cloud with a few simple clicks.
- > Easily automate the transfer of bioprocess data to the cloud. There is no need for manual uploads or exports using Excel sheets.

## Simple upload of historic data to the cloud

- > Easily upload historic bioprocess data from DASware control 5 or DASware control 6 software to the cloud.

## Upload of data from external analyzers

- > Upload data generated with external analyzers, for example cell density, viability, and nutrient concentrations.

## Upload of data from other SCADA software

- > Upload bioprocess data from BioCommand® software to the cloud.
- > Upload bioprocess data generated with non-Eppendorf bioprocess control software.

### Select finished runs in DASware control

Configuration

Plant

Workflows

Workflows

Templates

Edited

Waiting

Running

Finished

Inconsistent

Deleted

### Upload bioprocess data from finished runs in DASware control to the cloud

Plant Configuration

Workflows - Finished

Process Overview: LE493 DASbox glass 0.3L ExpiCHO batch pH optimum

Unit 1

Unit 2

Unit 4

Drag a column header and drop it here to group by that column

Started	Stopped	Title	Created	Comment
2024-05-28 09:33	2024-05-28 13:53	LE491 cleanup	2024-05-28 09:32	Cell Cultivation: pH Control using CO2 and Base
2024-05-22 10:09	2024-05-28 09:19	LE491 Pumpcalibration experiment 0.2L ExpiCHO ph optimum	2024-05-22 10:06	Cell Cultivation: pH Control using CO2 and Base
2024-05-21 10:39	2024-05-22 08:33	LE491_pH Kalibrierung und Sondentest	2024-05-21 10:32	Cell Cultivation: pH Control using CO2 and Base
2024-05-17 17:15	2024-05-17 17:17	dwcadministrator 7a52c7dd	2024-05-17 17:15	Fermentation of Aerobic Cells: DASbox: Control
2024-05-16 08:37	2024-05-17 12:32	0.3sc and 0.3c 1 SLPH profile to kLa	2024-05-16 08:37	Cell Cultivation: Controller not connected. Addit
2024-05-14 10:37	2024-05-15 12:02	Administrator b0937fda	2024-05-14 10:36	Cell Cultivation: pH Control using CO2 and Base
2024-05-08 16:22	2024-05-10 08:38	LE490 temperatur	2024-05-08 16:21	Cell Cultivation: pH Control using CO2 and Base
2024-05-07 10:50	2024-05-14 10:06	LE490_DASbox gla	2024-05-07 10:49	Cell Cultivation: pH Control using CO2 and Base
2024-05-05 15:47	2024-05-06 15:35	0.3sc and 0.3c 2 SLPH profile to kLa	2024-05-05 15:46	Cell Cultivation: Controller not connected. Addit
2024-05-03 14:53	2024-05-05 15:44	0.3sc and 0.3c 3 SLPH profile to kLa	2024-05-03 14:53	Cell Cultivation: Controller not connected. Addit
2024-04-30 11:36	2024-05-01 16:52	0.3sc and 0.3c 5 SLPH profile to kLa	2024-04-30 11:24	Cell Cultivation: Controller not connected. Addit
2024-04-23 16:14	2024-04-26 10:56	LE488_DASbox glass vs single-use 0.3L ExpiCHO batch demo with Fiona_BWS	2024-04-23 14:10	Cell Cultivation: pH Control using CO2 and Base
2024-03-07 12:00	2024-03-15 12:02	LE486_DASbox 0.3L ExpiCHO batch_GBMdemo_BWS	2024-03-07 11:56	Cell Cultivation: pH Control using CO2 and Base
2024-02-26 19:27	2024-03-01 10:37	LE484_DASbox 0.3L ExpiCHO batch_new vessel_pH_2_BWS	2024-02-26 19:22	Cell Cultivation: pH Control using CO2 and Base

# Build on your Data

Contextualize data across devices, runs, and sites - To optimize, transfer, scale-up

With BioNsight cloud you can easily synchronize and compare data from different batches, making it easier to identify patterns and anomalies across runs, accelerating your development process, and ensuring consistent quality and performance. Collaborate with colleagues, from researchers to production managers, regardless of their location. Discuss your process data when transferring a process to a different site, or compare performance between small-scale and bench-scale systems during scale-up.

## Align, compare, and analyze your bioprocess runs

BioNsight cloud lets you easily synchronize and compare all relevant bioprocess data tracks, both from internal sensors and external analyzers. It therefore helps identifying patterns and discovering interdependencies.

### Align data


Align all process records to inoculation time

### Select runs for analysis

Bioprocess runs of interest, easily selected from a drop-down list of all your historic runs

### Expand your data sets

Upload and display data generated with external analyzers, for example cell density, viability, and nutrient concentrations





# From Data to Insights

Transform bioprocessing with artificial intelligence

BioNsight cloud provides a seamless journey from data to insight through its integration with DataHowLab. DataHowLab is a unique bioprocess data analytics and modelling solution that empowers scientists to generate powerful AI-enabled process insights with their process data. Advanced technologies are embedded in a user-friendly environment developed for non-data-science experts.

## Key capabilities of DataHowLab include:

### Data analytics and visualizations

Effortlessly generate insights with DataHowLab's visualization boards. A wide variety of analyses and plots can be generated and customized in just a few clicks.

### Model-based process development

More than just data insight, models can be used for key development tasks such as process optimization, robustness analysis, optimal experimental design, and more.

### Access to advanced process models

Access a library of advanced process models, including DataHow's transformative AI-enabled hybrid models.

### Advanced analytics with ease

DataHowLab has been developed for process scientists, providing guidance with workflows and automations to make even the most advanced analysis possible.

Scientists are often unable to unlock the value of their data due access and analytical constraints. The partnership of BioNsight cloud and DataHowLab offers a transformative path to data-driven, digital process development.



## Native integration of BioNsight® cloud and DataHowLab

**Use data efficiently:** Easily send all your data from DASware control to DataHowLab via BioNsight cloud



**Save time:** Data from BioNsight cloud is pre-formatted. No need for manual processing for use with DataHowLab



**Gain convenience:** Single sign-on (SSO) to log in to DataHowLab with your BioNsight cloud user data



Find out more about DataHowLab  
[www.datahow.ch/products/datahowlab/](http://www.datahow.ch/products/datahowlab/)





## Find out more about BioNsight cloud

We'd love to tell you more about our software, answer any questions you may have, or schedule a software demo for you!

Visit our website for more information  
[www.eppendorf.group/BioNsight-cloud](http://www.eppendorf.group/BioNsight-cloud)



# Ordering Information

Description	Order No.
Software licence, for BioNsight® cloud access	7800001001
Software licence, for BioNsight® cloud access and DataHowLab connectivity	7800004001

**Your local distributor:** [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.group/BioNsight-cloud](http://www.eppendorf.group/BioNsight-cloud)

Microsoft® is a registered trademark of Microsoft Corp., USA. Eppendorf®, the Eppendorf Brand Design, and BioNsight® cloud are registered trademarks of Eppendorf SE, Germany. DASware® is a registered trademark of DASGIP Information and Process Technology GmbH, Germany. BioCommand® is a registered trademark of Eppendorf Inc., USA. All rights reserved, including graphics and images. Copyright © 2024 by Eppendorf SE.  
 GB3/PDF/0225/EBCC

