

# Enjoy Peace of Mind - Above and Beyond Sample Safety

## Get connected to VisioNize® Lab Suite

Whether you are looking for a remote monitoring solution, a lab management system or ways to easily access device performance data, VisioNize Lab Suite has all the services to help you build a smart laboratory with enhanced sample safety, productivity and reliability.



### For starters:

The Essential plan provides you with the basic functions to begin the digital management of your laboratory, starting with increasing your sample safety.

### For more peace of mind:

Create device documentation at the press of a button whilst benefiting from enhanced notification options. Enjoy the enriched feature set in VisioNize Lab Suite for more convenience!

	For Free (3 connected devices included*)	For more peace of mind (Billed annually per connected device**)
	VisioNize Lab Suite Essentials	VisioNize Lab Sample Safety+
<b>Real-Time Monitoring</b> Monitor connected lab devices and track status, parameters, and performance history from anywhere at any time.	✓	✓
<b>Instant Notifications</b> Receive alarm, error or event notifications to address any issues quickly.	Email only	Email & SMS messaging
<b>Non-Connectable Devices</b> Manage an unlimited number of non-connectable devices by manual data entry, e.g. serial number and location.	✓	✓
<b>Software Updates</b> Receive remote device software updates for connected Eppendorf lab devices.	✓	✓
<b>Task Management</b> Plan and assign device-related tasks and get notified when due.	✓	✓
<b>Document Management</b> Quickly retrieve relevant device documentation, e.g. operating manuals and other SOPs.	✓	✓
<b>Device Data Export</b> Export device parameter, event data and any user interactions.	✓	✓
<b>Lab Equipment Booking</b> Book and plan the usage of your devices in one shared calendar to boost productivity and smooth execution of your experiments, e.g. PCR run.	✓	✓
<b>Device Data Access</b> Quickly retrieve all relevant device data and go back in time with the history view.	7 days	90 days
<b>Audit Trail</b> Track user interactions or alarm acknowledgements with time-stamped data in an audit trail.	—	✓
<b>Response-Based Alarm Escalation</b> Define escalation schemes to individual or multiple recipient groups based on presence in the laboratory.	—	✓

\* Additional devices require device licenses which will be billed annually. Hardware and software upgrades are not included.

\*\* All pricing information are available online at: [www.eppendorf.com/visionize-subscription](http://www.eppendorf.com/visionize-subscription)



Contact us for more information: [www.eppendorf.com/visionize-subscription](http://www.eppendorf.com/visionize-subscription)

