

NEW: Conical Tube SnapTec® 50




Done in a Snap

Fast and safe sample handling with Eppendorf Conical Tubes SnapTec® 50

When things need to go fast, you feel like you need more than two hands. Not with the new Eppendorf Conical Tubes SnapTec® 50. They can be opened and closed quickly and safely with one hand – with either the left or the right – and enable significantly faster liquid withdrawal or addition. The cap is securely connected to the tube, eliminating the risk of mixing up lids while they lie on the bench.

Easily integrate this Conical Tube SnapTec 50 into your daily laboratory routine. With the dimensions comparable to a classic 50 mL conical tube, they fit perfectly into a broad range of your existing lab equipment.



Discover the SnapTec® 50 Cap Solution

Innovative SnapTec cap
Protected by two technical patents
(see also Eppendorf Conical Tubes 25 mL with SnapTec cap)

Nominal volume of 45 mL*

Autoclavable at 121°C, 20 min*

Centrifugation stability up to 17,000 x g*

Tightness of seal from -86°C up to 100°C*

The new Eppendorf Conical Tubes SnapTec® 50 are recommended for:

- > Centrifugation
- > Preparation of buffers, aliquots, working solutions
- > Sample preparation
- > Sample collection
- > Temperature incubation of samples

Single-handed operation
The SnapTec cap supports optimized and quicker handling of sample solutions, e.g. in multistep protocols

Reduced risk of contamination
The SnapTec cap avoids risk of mixing lids when working in multistep protocols

Implement immediately
Comparable dimension to 50 mL conical tubes with screw cap allows straightforward integration into a broad range of existing lab equipment*

Quality Certifications:

- guaranteed quality**
- sterile** (pyrogen-free, certified purity grade, sterile)
- PCR clean** (PCR inhibitor-free, certified purity grade, human DNA-free)

*Please note the instructions for use for details (www.eppendorf.com/manuals): Eppendorf Conical Tubes 25 mL and SnapTec® 50

Less risk

Reduced risk of contamination:
No need to put down the tube lid on the bench – the cap is securely connected to the tube – even when open. No risk of mixing up lids.

Single-handed operation

With the new SnapTec 50 conical tube reliably open and close the tube for quick liquid extraction or addition with only one hand.

Maximal comfort and safety

Learn more about the secure fitting cap, centrifugation stability, and leachable levels in our new application note*.
www.eppendorf.com/appnote445

*Eppendorf Conical Tubes SnapTec® 50: Fast and safe sample handling

Ordering information

Eppendorf Conical Tubes SnapTec® 50

Eppendorf Quality™, (240 tubes, 8x30 tubes/bag)

PCR clean, (240 tubes, 8x30 tubes/bag)

Sterile, (200 tubes, 10x20 tubes/bag)

Order no.

0030 118 650

0030 118 669

0030 118 677

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com/SnapTec50