

Assignment of epT.I.P.S.[®] Box 2.0 Options and Suitable Reload Variants

Box size		Color coding	Volumes	Reload variant stacked against each other	Reload variant stacked on top of each other	Sterile Reload
1		Dark grey	0.1 - 10 µL		■	
		Medium grey	0.1 - 20 µL	■		●
		Light grey	0.5 - 20 µL L		■	●
		Light pink	0.1 - 20 µL	■		
		Yellow	2 - 200 µL		■	●
		Light yellow	0.5 - 100 µL	■		
		Orange	20 - 300 µL	■		●
2		Blue	50 - 1,000 µL	■		●
		Green	50 - 1,250 µL	■		●
3		Dark green	50 - 1,250 µL L	■		●
		Red	0.25 - 2.5 mL	■		
		Violet	0.1 - 5 mL	Only in Racks, no Reloads available		

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · 22339 Hamburg · Germany eppendorf@eppendorf.com

www.eppendorf.com

Reload variants:

- = epT.I.P.S.[®] Standard
- = epT.I.P.S.[®] BioBased