

## Assignment of epT.I.P.S.<sup>®</sup> Reload and Rack Variants to Available Purity Grades

Reload variants:		epT.I.P.S. <sup>®</sup>			ep Dualfilter T.I.P.S. <sup>®</sup>				
■ = Reload									
● = Sterile Reload		<b>Eppendorf Quality</b>	<b>PCR clean</b>	<b>Biopur<sup>®</sup></b>	<b>PCR clean/Sterile</b>	<b>Forensic DNA Grade</b>	<b>SealMax<sup>®</sup> PCR clean/Sterile</b>	<b>SealMax<sup>®</sup> Biopur<sup>®</sup></b>	
▲ = Rack									
Color Coding	Volume								
Dark grey	0.1 – 10 µL	■	■		● (M <sup>1)</sup> ▲ (S, M <sup>1)</sup>	▲ (M <sup>1)</sup>	▲ (S)		
Medium grey	0.1 – 20 µL	■	■	● ▲					
Light grey	0.5 – 20 µL L	■	■		● ▲		▲	●	
Light pink	0.1 – 20 µL	■	■		▲				
Yellow	2 – 20 µL				● ▲	▲			
Yellow	2 – 100 µL				● ▲		▲	●	
Yellow	2 – 200 µL	■	■	● ▲	● ▲	▲	▲	●	
Light yellow	0.5 – 100 µL	■	■		▲				
Orange	20 – 300 µL	■	■	● ▲	● ▲		▲	●	
Blue	50– 1,000 µL	■	■	● ▲	● ▲	▲	▲	●	
Green	50 – 1,250 µL	■	■	● ▲					
Dark green	50 – 1,250 µL L	■	■	● ▲	● ▲				
Red	0.25 – 2.5 mL	■	■	▲	▲				
Violet	0.1 – 5 mL		▲		▲				
Violet	0.2 – 5 mL L		▲	▲	▲				
Turquoise	0.5 – 10 mL			▲					

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf SE · 22339 Hamburg · Germany [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)

1) Please note that 0.1-10 µL ep Dualfilter T.I.P.S.<sup>®</sup> size M are coloured in medium grey