



UK Health
Security
Agency

Certificate of Containment Testing

Containment Testing of 1500ml bottles in a R9A2 Rotor, Operated in as High-Speed Centrifuge

Report No. 23/013F

Report Prepared For: Eppendorf SE, Hamburg, Germany

Issue Date: 28 November 2023

Test Summary

The 1500ml bottle (part number 5721411035) in Rotor R9A2 (max. speed 8500 rpm) were containment tested at 8500 rpm in an Eppendorf Himac Technologies High-Speed centrifuge, using Annex AA of IEC 61010-2-020:2016 (3rd Ed.). The sealed bottles were shown to contain a spill.

Report Written By

Name: Ms Helen Hookway

Title: Biosafety Team Leader

Report Authorised By

Name: Mr Thomas Pottage

Title: Senior Biosafety Scientist

Please be aware that the use of the Royal Coat of Arms is highly restricted and cannot be copied. Please do not put the UKHSA logo on your website or use our name to endorse your products. Any reference to UKHSA needs to be approved by us before it can be used.