

Eppendorf Certificate

Certificate of Quality

Membrane Filter, 0.45µm for Easypet® 3
International Order No.: 4421 601.009


The above mentioned product, distributed by Eppendorf SE, is subject to careful quality inspection prior to release for distribution.

Each lot is tested 100% during manufacture for membrane and housing integrity by a leakage test under automatic conditions and by visual inspection.


Each lot is subject to further tests to assure that this product is only released, if it fulfills the requirements and meets the specifications of the technical data.

Sterility is confirmed according to ISO 11135.

Hamburg, July 2023



Monika Schneider
Vice President Global Quality Management
& Regulatory Affairs



Christian Behrens
Product Lifecycle Manager
Business Division Liquid Handling

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf®, the Eppendorf Brand Design and Easypet® are registered trademarks of Eppendorf SE, Germany.
All rights reserved, incl. graphics and images. Copyright ©2023 by Eppendorf SE.

www.eppendorf.com